

Product datasheet

Duo Advanced

An excellent 2-section treatment couch, with an emphasis on treatment in a seated position. Also particularly suitable as a robust base for general examination of flexion and extension functions.

Standard features

- Two-section couch with adjustable back section
- Back section supported by 2 gas springs for maximum stability
- Electric or hydraulic height adjustment
- Gymna Stability Profile (GSP)
- Cushions with 40 mm padding
- Breathing hole in back incl. cushion
- Central wheel-lifting mechanism with 4 casters
- Non slip feet
- Cushions in a choice of 11 colours
- Prepared for hip traction support

Ordering data

125 001	Duo Advanced peripheral switch
125 002	Duo Advanced foot switch
125 003	Duo Advanced hydraulic

Specifications

Length x width:	± 202 x 67cm
Height adjustment:	± 45 – 95 cm
Weight:	± 80 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Hydraulic pump:	10.000 N
Power supply:	100 – 240 VAC, 50-60 Hz
Current consumption:	1,5 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC
Guarantee:	3 years (upholstery 1 year)



Optional

121 404	Set brackets for fastening straps
121 402	Cushions 60 cm wide (Advanced)
121 403	Cushions 80 cm wide (Advanced)

Accessories

122 405	Paper roll holder, for undercarriage
125 651	Paper roll holder, for back section
121 409	Extra footswitch for height adjustment
125 652	Hand switch for height adjustment
337 018	Paper rolls (per 12 pieces)
107 819	Traction support for hip
121 401	Gynaecological leg support
333 245	Traction device AutoTrac 460
121 434	Attachment frame for AutoTrac 450/460
332 915	Foot end protector for cushion
xxx xxx	Terry-cloth cover, stretch material (available in 9 colours 'Style')





Gymna Stability Profile

Duo Advanced



Functionality/dimensions



Cushion colours

Style



Fashion

