

# Product datasheet

## Duo Luxe

*An excellent 2-section treatment table, with an emphasis on treatment in a seated position. Also particularly suitable as a robust base for general examination of flexion and extension functions.*

### Standard features

- Two-section table with adjustable back section.
- Gymna cushion heating with 2 temperature settings and hand control.
- Ergomax Face Support for optimum patient comfort while lying in prone position (one supplied as standard).
- Back section supported by 2 gas springs for maximum stability.
- Electrical height adjustment, operated by peripheral switch from all 4 sides of the table.
- Gymna Stability Profile (GSP).
- Cushions with duo-padding, with a high-density lower layer.
- The Luxe cushions are finished with tough saddle stitching.
- Central wheel-lifting mechanism with 4 casters.
- Non slip feet.
- Cushions in a choice of 11 colours.
- Prepared for hip traction support.

### Ordering data

125 004 Duo Luxe incl. peripheral switch and cushion heating

### Specifications

Length x width:	± 202 x 67cm
Height adjustment:	± 45 – 95 cm
Weight:	± 84 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Cushion temperature:	35 or 38 °C
Power supply:	100 – 240 VAC, 50-60 Hz
Current consumption:	3,0 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC
Guarantee:	3 years (upholstery 1 year)



### Optional

121 404	Set brackets for fastening straps
125 608	Cushions 60 cm wide
125 609	Cushions 80 cm wide

### Accessories

150 414	Ergomax Face Support
122 405	Paper roll holder, for undercarriage
125 651	Paper roll holder, for back section
337 018	Paper rolls (per 12 pieces)
121 409	Extra footswitch for height adjustment
107 819	Traction support for hip
121 401	Gynaecological leg support
333 245	Traction device AutoTrac 460
121 434	Attachment frame for AutoTrac 450/460
332 915	Foot end protector for cushion
xxx xxx	Terry-cloth cover, stretch material (available in 9 colours 'Style')



Gymna Stability Profile




Gymna cushion heating inside

**Duo Luxe**



**NEW:** Ergomax Face Support

**Functionality/dimensions**



**Cushion colours**

**Style**



**Fashion**

