

Product datasheet

Duoflex Advanced

The perfect solution for anyone who desires simplicity and quality above all. Perfect for use in massage and relaxation in lying positions. The armrests provide an extremely comfortable ventral position for your patient. Offering the best possible relaxation of the shoulder and neck muscles during your treatment.

Standard features

- Two-section couch with adjustable head section with arm supports
- Electric or hydraulic height adjustment
- Gymna Stability Profile (GSP)
- Cushions with 40 mm padding
- Breathing hole in head section, incl. cushion
- Central wheel-lifting mechanism with 4 casters
- Non slip feet
- Cushions in a choice of 11 colours
- Prepared for shoulder traction supports

Ordering data

125 031	Duoflex Advanced peripheral switch
125 032	Duoflex Advanced foot switch
125 033	Duoflex Advanced hydraulic

Specifications

Length x width:	± 202 x 67cm
Height adjustment:	± 45 – 95 cm
Weight:	± 80 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Hydraulic pump:	10.000 N
Power supply:	100 – 240 VAC, 50-60 Hz
Current consumption:	1,5 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC
Guarantee:	3 years (upholstery 1 year)



Optional

125 610	Motorized head section
121 404	Set brackets for fastening straps
121 402	Cushions 60 cm wide (Advanced)
121 403	Cushions 80 cm wide (Advanced)

Accessories

122 405	Paper roll holder, for undercarriage
337 018	Paper rolls (per 12 pieces)
121 409	Extra footswitch for height adjustment
125 652	Hand switch for height adjustment
107 819	Traction support for shoulders
121 401	Gynaecological leg support
333 245	Traction device AutoTrac 460
121 434	Attachment frame for AutoTrac 450/460
332 915	Foot end protector for cushion
xxx xxx	Terry-cloth cover, stretch material (available in 9 colours 'Style')





Gymna Stability Profile

Duoflex Advanced



Functionality/dimensions



Cushion colours

Style



Fashion

