

# Product datasheet

## Duoflex Luxe

*The perfect solution for anyone who desires simplicity and quality above all. Perfect for use in massage and relaxation in lying positions. The armrests provide an extremely comfortable ventral position for your patient. Offering the best possible relaxation of the shoulder and neck muscles during your treatment.*

### Standard features

- Two-section table with adjustable head section with arm supports.
- Gymna cushion heating with 2 temperature settings and hand control.
- Ergomax Face Support for optimum patient comfort while lying in prone position (one supplied as standard).
- Electrical height adjustment, operated by peripheral switch from all 4 sides of the table.
- Gymna Stability Profile (GSP).
- Cushions with duo-padding, with a high-density lower layer.
- The Luxe cushions are finished with tough saddle stitching.
- Central wheel-lifting mechanism with 4 casters.
- Non slip feet.
- Cushions in a choice of 11 colours.
- Prepared for shoulder traction supports.

### Ordering data

125 034 Duoflex Luxe incl. peripheral switch and cushion heating

### Specifications

Length x width:	± 202 x 67cm
Height adjustment:	± 45 – 95 cm
Weight:	± 84 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Cushion temperature:	35 or 38 °C
Power supply:	100 – 240 VAC, 50-60 Hz
Current consumption:	3,0 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC
Guarantee:	3 years (upholstery 1 year)



### Optional

125 610	Motorized head section
121 404	Set brackets for fastening straps
125 608	Cushions 60 cm wide
125 609	Cushions 80 cm wide

### Accessories

150 414	Ergomax Face Support
122 405	Paper roll holder, for undercarriage
337 018	Paper rolls (per 12 pieces)
107 819	Traction support for shoulders
121 401	Gynaecological leg support
333 245	Traction device AutoTrac 460
121 434	Attachment frame for AutoTrac 450/460
332 915	Foot end protector for cushion
xxx xxx	Terry-cloth cover, stretch material (available in 9 colours 'Style')



Gymna Stability Profile




Gymna cushion heating inside

## Duoflex Luxe



**NEW:** Ergomax Face Support

### Functionality/dimensions



### Cushion colours

#### Style



#### Fashion

