

# Product datasheet

## Duoplan Luxe

*An excellent 2-section treatment table with individually adjustable leg sections. Emphasis on treatment in a seated position and Particularly suitable for treatment of the lower limbs, for example for lymph drainage, traumatology or podotherapy*

### Standard features

- Two-section table with adjustable back section and split leg section.
- Gymna cushion heating with 2 temperature settings and hand control.
- Ergomax Face Support for optimum patient comfort while lying in prone position (one supplied as standard).
- All sections supported by 2 gas springs.
- Electrical height adjustment, operated by peripheral switch from all 4 sides of the table.
- Gymna Stability Profile (GSP).
- Cushions with duo-padding, with a high-density lower layer.
- The Luxe cushions are finished with tough saddle stitching.
- Central wheel-lifting mechanism with 4 casters.
- Non slip feet.
- Cushions in a choice of 11 colours.
- Prepared for hip traction support.

### Ordering data

125 014 Duoplan Luxe incl. peripheral switch and cushion heating

### Specifications

Length x width:	± 202 x 67cm
Height adjustment:	± 45 – 95 cm
Weight:	± 95 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Cushion temperature:	35 or 38 °C
Power supply:	100 – 240 VAC, 50-60 Hz
Current consumption:	3,0 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC
Guarantee:	3 years (upholstery 1 year)



### Optional

- 121 406 Set brackets for fastening straps
- 125 608 Cushions 60 cm wide

### Accessories

- 150 414 Ergomax Face Support
- 122 405 Paper roll holder, for undercarriage
- 125 651 Paper roll holder, for back section
- 337 018 Paper rolls (per 12 pieces)
- 107 819 Traction support for hip
- 332 915 Foot end protector for cushion
- xxx xxx Terry-cloth cover, stretch material (available in 9 colours 'Style')



Gymna Stability Profile




Gymna cushion heating inside

**Duoplan Luxe**



**NEW:** Ergomax Face Support

**Functionality/dimensions**



Option 60 cm

**Cushion colours**

**Style**

-   
*Pacific*
-   
*Carbon*
-   
*Chrome*
-   
*Ivory*
-   
*Sahara*
-   
*Marble*
-   
*Havana*
-   
*Terracotta*

**Fashion**

-   
*Agave*
-   
*Mandarin*
-   
*Atoll*