

Product datasheet

Duoplus Luxe

The perfect solution for anyone who desires simplicity and quality above all. Perfect for use in massage and relaxation in lying positions

Standard features

- Two-section table with adjustable head section (gas spring supported).
- Gymna cushion heating with 2 temperature settings and hand control.
- Ergomax Face Support for optimum patient comfort while lying in prone position (one supplied as standard).
- Electrical height adjustment, operated by peripheral foot switch from all 4 sides of the table.
- Gymna Stability Profile (GSP).
- Cushions with duo-padding, with a high-density lower layer.
- The Luxe cushions are finished with tough saddle stitching.
- Central wheel-lifting mechanism with 4 casters.
- Non slip feet.
- Cushions in a choice of 11 colours.
- Prepared for shoulder traction supports.

Ordering data

125 024 Duoplus Luxe incl. peripheral switch and cushion heating

Specifications

| | |
|------------------------|--|
| Length x width: | ± 202 x 67cm |
| Height adjustment: | ± 45 – 95 cm |
| Weight: | ± 82 kg |
| Max. lifting capacity: | 200 kg |
| Electric motor: | 6.000 N |
| Cushion temperature: | 35 or 38 °C |
| Power supply: | 100 – 240 VAC, 50-60 Hz |
| Current consumption: | 3,0 A |
| Insulation class: | II |
| Moister resistance: | IP 21 |
| Classification: | I |
| Conformity: | IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC |
| Guarantee: | 3 years (upholstery 1 year) |



Optional

| | |
|---------|-----------------------------------|
| 125 610 | Motorized head section |
| 121 404 | Set brackets for fastening straps |
| 125 608 | Cushions 60 cm wide |
| 125 609 | Cushions 80 cm wide |

Accessories

| | |
|---------|--|
| 150 414 | Ergomax Face Support |
| 122 405 | Paper roll holder, for undercarriage |
| 337 018 | Paper rolls (per 12 pieces) |
| 107 819 | Traction support for shoulders |
| 121 401 | Gynaecological leg support |
| 333 245 | Traction device AutoTrac 460 |
| 121 434 | Attachment frame for AutoTrac 450/460 |
| 332 915 | Foot end protector for cushion |
| xxx xxx | Terry-cloth cover, stretch material (available in 9 colours 'Style') |



Gymna Stability Profile





Gymna cushion heating inside

Duoplus Luxe



NEW: Ergomax Face Support

Functionality/dimensions



Cushion colours

Style



Fashion

