

# Product datasheet

## Osteoflex Luxe

Particularly suitable for those who concentrate on treatment of the spinal column. Thanks to the unique fold-down side supports, you can stand closer to the patient. This enables you to exert more force with less strain on your own shoulders and back. The arm supports provide a very comfortable prone position, with the best possible relaxation of the shoulder and neck muscles. Perfectly suited for osteopathic manipulations.

### Standard features

- Three-section table with adjustable head-middle- and leg section.
- Gymna cushion heating with 2 temperature settings and hand control.
- Ergomax Face Support for optimum patient comfort while lying in prone position (one supplied as standard).
- All sections gas spring supported.
- Arm supports in the head section.
- Unique fold-down supports in the middle section.
- Electrical height adjustment, operated by peripheral foot switch from all 4 sides of the table.
- Gymna Stability Profile (GSP).
- Cushions with duo-padding, with a high-density lower layer.
- The Luxe cushions are finished with tough saddle stitching.
- Central wheel-lifting mechanism with 4 casters.
- Non slip feet.
- Cushions in a choice of 11 colours.
- Prepared for hip traction support.
- Prepared for shoulder traction supports.

### Ordering data

125 064 Osteoflex Luxe incl. peripheral switch and cushion heating

### Specifications

Length x width:	± 202 x 67cm
Height adjustment:	± 45 – 95 cm
Weight:	± 110 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Cushion temperature:	35 or 38 °C
Power supply:	100 – 240 VAC, 50-60 Hz
Current consumption:	3,0 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC
Guarantee:	3 years (upholstery 1 year)



### Optional

125 602	Electrical adjustable middle section
125 610	Motorized head section (not in combination with 125 602)
121 405	Set brackets for fastening straps
125 608	Cushions 60 cm wide (Luxe)

### Accessories

150 414	Ergomax Face Support
122 405	Paper roll holder, for undercarriage
337 018	Paper rolls (per 12 pieces)
107 819	Traction support for hips or shoulders
121 401	Gynaecological leg support
333 245	Traction device AutoTrac 460
121 434	Attachment frame for AutoTrac 450/460
332 915	Foot end protector for cushion
xxx xxx	Terry-cloth cover, stretch material (available in 9 colours 'Style')



Gymna Stability Profile





Gymna cushion heating inside

**Osteoflex Luxe**



**NEW:** Ergomax Face Support

**Functionality/dimensions**



**Cushion colours**

**Style**



**Fashion**

