

Product datasheet

Trioplus Luxe

A 3-section table in the purest form – stable and functional. Adjustable headrest, middle- and leg section. Extremely versatile and suited for treatment in sitting, lying, lumbar flexion and Trendelenburg position. Guarantees years of problem-free support for your manipulations.

Standard features

- Three-section table with adjustable head-, middle-, and leg section.
- Gymna cushion heating with 2 temperature settings with hand control.
- Ergomax Face Support for optimum patient comfort while lying in prone position (one supplied as standard).
- All sections gas spring supported.
- Electrical height adjustment, operated by peripheral foot switch from all 4 sides of the table.
- Gymna Stability Profile (GSP).
- Cushions with duo-padding, with a high-density lower layer.
- The Luxe cushions are finished with tough saddle stitching.
- Central wheel-lifting mechanism with 4 casters.
- Non slip feet.
- Cushions in a choice of 11 colours.
- Prepared for hip traction support.
- Prepared for shoulder traction supports.

Ordering data

125 044 Trioplus Luxe incl. peripheral switch and cushion heating

Specifications

Length x width:	± 202 x 67cm
Height adjustment:	± 45 – 95 cm
Weight:	± 101 kg
Max. lifting capacity:	200 kg
Electric motor:	6.000 N
Cushion temperature:	35 or 38 °C
Power supply:	100 – 240 VAC, 50-60 Hz
Current consumption:	3,0 A
Insulation class:	II
Moister resistance:	IP 21
Classification:	I
Conformity:	IEC/EN 60601-1 IEC/EN 60601-2 CAN/CSA C22.2 No. 601.1 MDD 93/42/EEC
Guarantee:	3 years (upholstery 1 year)



Optional

125 602	Electrical adjustable middle section
125 610	Motorized head section (not in combination with 125 602)
121 404	Set brackets for fastening straps
125 608	Cushions 60 cm wide
125 609	Cushions 80 cm wide

Accessories

150 414	Ergomax Face Support
122 405	Paper roll holder, for undercarriage
337 018	Paper rolls (per 12 pieces)
107 819	Traction support for hips or shoulders
121 401	Gynaecological leg support
333 245	Traction device AutoTrac 460
121 434	Attachment frame for AutoTrac 450/460
332 915	Foot end protector for cushion
xxx xxx	Terry-cloth cover, stretch material (available in 9 colours 'Style')



Gymna Stability Profile





Gymna cushion heating inside

Trioplus Luxe



NEW: Ergomax Face Support

Functionality/dimensions



Cushion colours

Style



Fashion

